

First Aid Kit Contents:

Protective gloves	2 pair
1" x 3" Fabric Bandage	50 ea.
Knuckle fabric bandage	40 ea.
Antiseptic spray (aerosol)	1 ea.
Eye and skin flushing solution 8 oz	1 ea.
Eye wash 4 oz.	1 ea.
Burn gel	20 ea.
Roller gauze 1"	1 ea.
Roller gauze 2"	1 ea.
Roller gauze 3"	1 ea.
2"x2" sterile gauze pad	6 ea.
3" x 3" sterile gauze pad	4 ea.
4" x 4" sterile gauze pad	4 ea.
Gauze compress 36" x 36"	1 ea.
Eye pad w/adhesive strip	4 ea.
Cold pack	1 ea.
Adhesive tape 3 cut	1 ea.
Tri-bandage	1 ea.
Microshield clear mouth barrier	1 ea.
Scissors	1 pr.
Tweezers	1 pr.
First Aid manual	1 ea.

Automated External Defibrillator Maintenance Checklist

Date _____ Location _____
 Inspection Performed by _____

Criteria	Status											
	January	February	March	April	May	June	July	August	September	October	November	December
AED												
Visible placement												
Check expiration dates of pads/battery in AED, replace if necessary.												
Check the status indicator light												
Note absence of visual/audible service alarm												
Check unit and accessories for damage, dirt, and contamination. Clean or replace as necessary.												
Date and Initial AED Inspection Tag												
SUPPLIES												
Two sets of AED pads in sealed package												
Check expiration date on extra pads package												
Check expiration date on extra battery package												
CPR shield with one-way valve												
Examination gloves												
Razor/Scissor												

■ Please refer to manufacturer's User's Manual for more information and proper annual maintenance procedures.